

ADS EXPRESS

PORTABLE ANESTHESIA MACHINE



Rugged, compact and light weight, the ADS Express is ready for the demanding emergency surgery environment. Weighing in at less than 12 Kilos the ADS express is ready to go anywhere. The ADS Express features hypoxia-safe N₂O and O₂ flowmeters, APL valve, O₂ Flush and built-in quick fill 1000ml CO₂ absorber.

The ADS Express is well suited for:

- Office Based Anesthesia
- Remote Field Hospitals
- Disaster Recovery
- Dental Anesthesia
- Military Field Hospitals
- Mobile - Aircraft / Ships / Ambulance



Advanced, Affordable Medical Technology For The World

ADS EXPRESS - PORTABLE ANESTHESIA MACHINE TECHNICAL SPECIFICATIONS

Patient scope: Pediatric (3kg) to Adult
Tidal Volume: 20-1600 ML

Working pressure:

O₂ 0.3-0.4 Mpa
N₂O 0.3-0.4 Mpa

Flowmeter Range:

O₂ 0.1-10 LPM
N₂O 0.1-10 LPM
ORC Hypoxic Guard O₂ > 25%

Flowmeter Range (continue):

APL valve range: 0.5 – 6.0 Kpa

O₂ Flush: 30-60 L/min

NIST type O₂ and N₂O fittings

Selectatec type vaporizer:

Compatible with Sevoflurane,
Enflurane, Isoflurane, Halothane


Optional:

InfiVent Ventilator

Mobile Cart



3 Side Ports
(Expiratory, Inspiratory & Exhaust)



2 Back Side Inlets
(O₂ & N₂O)

CE

Authorized Distributor:



12151 62nd St N #5 • Largo FL 33773 • USA

Phone: (1) 727-531-8434 • Fax: (1) 727-531-8436

Web: www.infiniummedical.com • Email: sales@infiniummedical.com